GE Healthcare

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## 参考文献 (Reference)

1. Cardiac screening examination of the fetus: guidelines for performing the 'basic' and 'extended basic' cardiac scan. *Ultrasound Obstet Gynecol* 2006; 27: 107-113
2. Haak MC, Twisk JW, VAN VUGT JM. How successful is fetal echocardiographic examination in the first trimester of pregnancy? *Ultrasound Obstet Gynecol* 2002; 20: 9-13.
3. Gregory R. DeVore, M.D. , et al. High Frequency Endovaginal Fetal Echocardiography. GE Healthcare.
4. Vinals F, et al. Fetal echocardiography at 11+0 to 13+6 weeks using four-dimensional spatiotemporal image correlation telemedicine via an Internet link: a pilot study. *Ultrasound Obstet Gynecol* 2008; 31: 633-638.
5. Sklansky MS, Nelson TR, Pretorius DH. Three-dimensional fetal echocardiography: gated versus nongated techniques. *J Ultrasound Med* 1998 Jul; 17(7): 451-76.
6. Shih J. C. ,et al.'Big-eyed frog'sign on spatiotemporal image correlation(STIC) in the antenatal diagnosis of transposition of the great arteries. *Ultrasound Obstet Gynecol* 2008; 32: 762-768.
7. American Institute of Ultrasound in Medicine. AIUM practice guideline for the performance of an antepartum obstetric ultrasound examination. *J Ultrasound Med* 2003; 22: 1116-1125.
8. American College of Obstetricians and Gynecologists. ACOG practice bulletin 58: ultrasonography in pregnancy. *Ultrasound Obstet Gynecol* 2004; 104: 1449-1458.
9. Gregory R. DeVore, MD, Bardo Polanko, BS. Tomographic Ultrasound Imaging of the Fetal Heart. *J Ultrasound Med* 2005; 24: 1685-1696.
10. Paladini D, et al. The role of spatio-temporal image correlation (STIC) with a tomographic ultrasound imaging (TUI) in the sequential analysis of fetal congenital heart disease. *Ultrasound Obstet Gynecol* 2006; 27: 555-561.
11. Turan S, et al. Standardization of the first-trimester fetal cardiac examination using spatiotemporal image correlation with tomographic ultrasound and color Doppler imaging. *Ultrasound Obstet Gynecol* 2009; 33: 652-656.

12. Lami Y, Roberto R, et al. STAR (Simple Targeted Arterial Rendering) Technique: a Novel and Simple Method to Visualize the Fetal Cardiac Outflow Tracts. *Ultrasound Obstet Gynecol* 2010 Sep 27.
13. Lami Y, Roberto R, et al. FAST (Four chamber view And Swing Technique) Echo: a Novel and Simple Algorithm to Visualize Standard Fetal Echocardiographic Planes. *Ultrasound Obstet Gynecol* 2010 Sep 27.
14. Yagel S, et al. Rendering in fetal cardiac scanning: the intracardiac septa and the coronal atrioventricular valve planes. *Ultrasound Obstet Gynecol* 2006; 28: 266-274.
15. Goncalves LF, Espinoza J, Lee W, et al. A new approach to fetal echocardiography: digital casts of the fetal cardiac chambers and great vessels for detection of congenital heart disease. *J Ultrasound med.* 2005; 24: 415-24.
16. Falkensammer P. Spatio-temporal imaging correlation for volume ultrasound- Study of the fetal heart. GE Healthcare: Zipf, Austria, 2005.
17. Pooh RK. New application of B-flow sonoangiography in perinatology. *Ultrasound Obstet Gynecol* 2000; 15: 163.
18. Voipe P, Campobasso G, Stanziano A, et al. Novel application of 4D sonography with B-flow imaging and spatiotemporal image correlation (STIC) in the assessment of the anatomy of pulmonary arteries in fetuses with pulmonary atresia and ventricular septal defect. *Ultrasound Obstet Gynecol* 2006; 28: 40-46.
19. Voipe P, Campobasso G, De Robertis V, et al. Two and four dimensional echo-cardiography with B-flow imaging and spatiotemporal image correlation in prenatal diagnosis of isolated total anomalous pulmonary venous connection. *Ultrasound Obstet Gynecol* 200; 30: 830-837.
20. Volpe P, et al. Fetal interrupted aortic arch: 2D-4D echocardiography associations and outcome. *Ultrasound Obstet Gynecol* 2010; 35: 302-309.
21. Yagel S, et al. The fetal venous system, Part II: ultrasound evaluation of the fetus with congenital venous system malformation or developing circulatory compromise. *Ultrasound Obstet Gynecol* 2010; 36: 93-111
22. Messing B, et al. Fetal cardiac ventricle volumetry in the second half of gestation assessed by 4D ultrasound using STIC combined with inversion mode. *Ultrasound Obstet Gynecol* 2007; 30: 142-151.
23. Rusk Ra, Mori Y, Davies CH, Irvine T, Kenny A, Sahn DJ. Comparison of ventricular volume and mass measurements from B- and C-scan images with the use of real-time 3-dimensional echocardiography: studies in an in vitro model. *J Am Soc Echocardiogr* 2000; 13: 910-917.

24. Schneider C, et al. Development of Z-scores for fetal cardiac dimensions from echocardiography Ultrasound Obstet Gynecol 2005; 26: 599-605.
25. Beryl R. Benacerraf, Carol B. Benson, Alfred Z. Abuhamad et al. Three and 4-Dimensional Ultrasound in obstetrics and Gynecology, Processing of the American Institute of Ultrasound in Medicine Consensus Conference. J Ultrasound Med 2005; 24: 1587-1597.
26. Luís F. Gonçalves et al, A Systematic Approach to Prenatal Diagnosis of Transposition of the Great Arteries Using 4-Dimensional Ultrasonography With Spatiotemporal Image Correlation. J Ultrasound Med 2004; 23:1225-1231
27. Yagel S, Shen D The added value of 3D/4D ultrasound imaging in fetal cardiology: has the promise been fulfilled Ultrasound Obstet Gynecol 2010; 35: 260-262
28. Garne E, et al. Prenatal diagnosis of severe structural congenital malformations in Europe. Ultrasound Obstet Gynecol 2005; 25: 6-11
29. Cardiac screening examination of the fetus: guidelines for performing the 'basic' and 'extended basic' cardiac scan. Ultrasound Obstet Gynecol 2006; 27: 107-113
30. Lee W, et al. ISUOG consensus statement: what constitutes a fetal echocardiogram? Ultrasound Obstet Gynecol 2008; 32: 239-242
31. Lee W, et al. The 'starfish' sign: a novel sonographic finding with B-flow imaging and spatiotemporal image correlation in a fetus with total anomalous pulmonary venous return. Ultrasound Obstet Gynecol 2010; 35: 124-125
32. B. TUTSCHEK and D. J. SAHN. Semi-automatic segmentation of fetal cardiac cavities: progress towards an automated fetal echocardiogram. Ultrasound Obstet Gynecol 2008; 32: 176-180
33. Jimmy Espinoza, et al. Prenatal Diagnosis of Coarctation of the Aorta With the Multiplanar Display and B-Flow Imaging Using 4-Dimensional Sonography. J Ultrasound Med 2009; 28: 1375-1378.
34. Yagel S, et al. 3D and 4D ultrasound in fetal cardiac scanning: a new look at the fetal heart. Ultrasound Obstet Gynecol 2007; 29: 81-95

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